

Product datasheet

Pulson 100

The Pulson 100, a basic ultrasound device

The Gymna Pulson 100 is a user-friendly, portable device with all basic functions. Equipped with multifrequency ultrasound heads and push buttons control for all functions.



Article number: 320.200

Characteristics

- Ultrasound :
 - Multifrequency head (1 and 3 MHz), 4 cm²
 - Detachable head
 - 1 ultrasound output connector
 - Pulsed and continuous mode (10 – 20 – 30 – 40 – 50 - 100%)
 - Acoustic and visual contact control
- Functionalities :
 - Indications: 23
 - Free mode: 1

Technical specifications

Mains voltages	: 100- 240-VAC, 50/60 Hz +/- 10 %
Maximum power -in operation	: 29 VA
Dimensions (b x h x d)	: 268 x 124 x 295 mm
Weight	: ca 3,150 kg
Safety	Class I
Isolation	: Type BF
MDD classification	: IIa
Safety examination	: CE from conformity with Directives MDD 93/42/EEC

Standard accessories

100 689	Mains lead
336 589	US-head, 1 & 3 MHz, ERA 4 cm ²
341 088	Contact gel, 500 ml
117 277	Sets of photos US-head placement
336 787	Indications leaflet

Manuals

323 011	Safety Instructions
331 958	User Manual English
376 126	CD-ROM User Manuals multi language

Optional accessories

336 578	US-head, 1 & 3 MHz, ERA 1 cm ²
302 955	Carrying bag 200 series
341 099	Contact-gel, 5L
341 121	Pump, 5L
320 804	Gymna Mobile

Therapies	
Ultrasound therapy (1&3 MHz)	■
User-interface	
Customized LCD for all parameters	■
Push Buttons for the important functions	■
Indications leaflet	■
Ultrasound therapy	
Hybrid treatment head 4 cm ² , 1&3 MHz	■
Hybrid treatment head 1 cm ² , 1&3 MHz	○
Acoustic and visual contact control led	■



- = Standard
- = Optional