


SPECIFICATIONS

| | | | |
|----------------------------|---------------------------|-----------------------------------------------------------------------------------|------------------------------------------------------------------------------|
| Dimensions: | 206 x 67 cm | Motor force: | 8 kN |
| Weight: | 104 kg | Motor usage: | Max. 2 minutes/ 18 minutes off |
| Height adjustment: | ± 47 – 97 cm | Hydraulic pump: | 10 kN |
| Lifting capacity: | Max. 200 kg | Frame material: | Steel frames |
| Mains voltage: | 100 - 240 VAC, 50 - 60 Hz | Frame finish: | Powder coated white RAL 9016 |
| Current: | Max. 3.0 A | Conformity: | IEC/EN 60601-1 |
| Standby power consumption: | Max. 3 W |  | IEC/EN 60601-1-2 |
| Insulation: | Class II | | EC directives 2017/745 (MDR) |
| Device Classification | I | | 2012/19/EC |
| Applied parts | Type BF | | 2011/65/EU |
| Humidity proof: | IP21 | Warranty: | Table 3 years, Premium upholstery 3 years, Classique upholstery 1 year |

Note: motor adjustment is slower when running on battery.

CUSHIONS UPHOLSTERY

| | |
|------------|-----------------------------------------------------------------------------|
| Classique: | PVC, phthalate free, biocompatible, abrasion resistant and flame retardant. |
| Premium: | PU, phthalate free, biocompatible, abrasion resistant and flame retardant. |

CUSHIONS FINISH

| | |
|--------------|-------------------------------------------------------------|
| Linear: | Single padded, 40 kg/m ³ |
| Progressive: | Dual padded, 100 kg/m ³ and 40 kg/m ³ |

OPTIONS

| | | | |
|-------------------|---------------------------------------------------------------------------------------|------------------------|---------------------------------------------------------------------------------------------------------------------------|
| Ergomax: | PU foam, closed skin, low density 5° Shore | Memory function: | Three programmable settings for height adjustment |
| Cushion heating: | 34 °C, 36 °C and 38 °C Bluetooth wireless operation | LED Pack: | Bluetooth wireless operation Chassis illumination 14 colours, including brightness control Power consumption = 10 W |
| Battery function: | Automatic timer 30 minutes Maximum consumption = 180 W Li-ion 25,9 VDC, 2,25 Ah | Safe Motion Detection: | Power consumption = 1 W |

FUNCTIONALITY/ DIMENSIONS

