

## SPECIFICATIONS

|                            |                           |                 |                                                                                    |
|----------------------------|---------------------------|-----------------|------------------------------------------------------------------------------------|
| Dimensions:                | 206 x 67 cm               | Motor force:    | 8 kN                                                                               |
| Weight:                    | 87 kg                     | Motor usage:    | Max. 2 minutes/ 18 minutes off                                                     |
| Height adjustment:         | ± 47 – 97 cm              | Hydraulic pump: | 10 kN                                                                              |
| Lifting capacity:          | Max. 200 kg               | Frame material: | Steel frames                                                                       |
| Mains voltage:             | 100 - 240 VAC, 50 - 60 Hz | Frame finish:   | Powder coated white RAL 9016                                                       |
| Current:                   | Max. 3.0 A                | Conformity:     | IEC/EN 60601-1                                                                     |
| Standby power consumption: | Max. 3 W                  |                 | IEC/EN 60601-1-2                                                                   |
| Insulation:                | Class II                  |                 | EU Regulation 2017/745                                                             |
| Device Classification      | I                         |                 | 2012/19/EU WEEE                                                                    |
| Applied parts              | Type BF                   |                 | 2011/65/EU RoHS                                                                    |
| Humidity proof:            | IP21                      | Warranty:       | Table 3 years,<br>Premium upholstery 3 years, Classique<br>2025 upholstery 2 years |



Note: motor adjustment is slower when running on battery.

## CUSHIONS UPHOLSTERY

|            |                                                                             |
|------------|-----------------------------------------------------------------------------|
| Classique: | PVC, phthalate free, biocompatible, abrasion resistant and flame retardant. |
| Premium:   | PU, phthalate free, biocompatible, abrasion resistant and flame retardant.  |

## CUSHIONS FINISH

|              |                                                             |
|--------------|-------------------------------------------------------------|
| Linear:      | Single padded, 40 kg/m <sup>3</sup>                         |
| Progressive: | Dual padded, 100 kg/m <sup>3</sup> and 40 kg/m <sup>3</sup> |

## OPTIONS

|                  |                                                           |                        |                                                                                           |
|------------------|-----------------------------------------------------------|------------------------|-------------------------------------------------------------------------------------------|
| Ergomax:         | PU foam, closed skin, low density 5° Shore                | Memory function:       | Three programmable settings for height adjustment                                         |
| Cushion heating: | 34 °C, 36 °C and 38 °C<br>Bluetooth wireless operation    |                        | Bluetooth wireless operation                                                              |
|                  | Automatic timer 30 minutes<br>Maximum consumption = 180 W | LED Pack:              | Chassis illumination 14 colours, including brightness control<br>Power consumption = 10 W |
|                  |                                                           | Safe Motion Detection: | Power consumption = 1 W                                                                   |

FUNKTIONALITÄT / ABMESSUNGEN

