

Product datasheet

Gymna Care 300

Powered by Fisiowarm

The Gymna Care 300 TECAR therapy device

The Gymna Care 300 is an advanced and user-friendly TECAR therapy device.

Effective and complete TECAR therapy to adapt perfectly to the specific patient needs.



Article number: 349 008

Characteristics

TECAR :

- 5 Frequencies
300kHz – 400kHz – 450kHz – 500kHz – 1 MHz
- AUTO - automatic scanning of pre-set frequencies

Emissions:

- Continuous
- Sensitive
- A-thermia 1/2/3

Operating modes:

- Unipolar
- Bipolar
- Neutro Dinamico
- Static Automatic

Energy transfer:

- Capacitive
- Resistive

Technical specifications Gymna Care 300 device:

Mains voltages: 110 – 230 VAC, 50 – 60Hz

Nominal power: 400 Watt

Dimensions: 288 (b) x 91 (h) x 308 (d) mm

Weight:

Gymna Care 300 4,65 kg, incl. iPad – Trolley 31,5 kg

Safety: Class I

Isolation: Type BF

MDD classification: IIa

User interface

- Gymna Care 300 iPad & App
- OLED display and rotary encoder with multi-function button

Standard packaging:

- Gymna Care 300 device
- iPad
- 349 107 - Unipolar kit (1 hand piece, 3 capacitive and 3 resistive electrodes)
- Neutral plate
- User manual
- Conductive cream (150 ml)
- 1 Kit of choice*
*Choose between the following kits:
- 349 118 - Bipolar kit (1 hand piece, 3 capacitive and 1 resistive electrode)
- Neutro Dinamico kit
- 349 129 Option 1 – cable and ND-handpiece for unipolar electrodes and 3 resistive unipolar electrodes
- 349 140 Option 2 – cable and handheld ND electrode
- 349 151 - Static automatic (cable, main cable and connections, 4 resistive clips electrodes, 2 rubber belts 130 cm)

Optional Accessories:

- 349 162 Cable for disposable neutral electrode
- 349 272 Disposable neutral electrodes (25 pcs)
- 349 173 Endoprobe unipolar
- 349 184 Endoprobe bipolar
- 349 195 Unipointer
- 349 206 Fascial full electrode with pouch
- 349 217 Facelift bipolar resistive electrode
- 349 228 Carrying bag (with wheels)
- 349 239 Trolley
- 349 250 Conductive cream 1000 ml