


SPECIFICATIONS

| | | | |
|------------------------|--|---|--|
| Dimensions : | 205 x 67 cm | Motor force: | 6 kN |
| Weight : | ±79 kg | Motor usage: | Max. 2 minutes on / 18 minutes off |
| Height adjustment : | 50 – 100 cm electric 52 – 99 cm hydraulic | Hydraulic pump: | 10kN |
| Lifting capacity: | Max. 200kg | Frame material: | Steel frames |
| Mains voltage: | 230 VAC, 50/60 Hz 120 VAC, 60 Hz 100 VAC, 50/60 Hz | Frame finish: | Powder coated white RAL 9016 |
| Power consumption: | Max. 2.0 A | Conformity: | IEC/EN 60601-1 IEC/EN 60601-1-2 EC directives 2017/745 (MDR) 2012/19/EC 2011/65/EU |
| Insulation: | Class II |  | Warranty: Table 3 years, Classique upholstery, 1 year |
| Device Classification: | I | | |
| Applied parts: | Type BF | | |
| Humidity proof: | IP21 | | |

CUSHION UPHOLSTERY

Classique: PVC, phthalate free, biocompatible, abrasion resistant and flame retardant.

CUSHIONS FINISH

Linear: Single padded, 40 kg/m³

OPTIONS

Wheels
Safety switch for foot- or hand control

Face opening
Ergomax face opening

DIMENSIONS / FUNCTIONALITY

