

## SPECIFICATIONS

Dimensions :	205 x 67 cm
Weight :	±81 kg
Height adjustment :	50 – 100 cm electric 52 – 99 cm hydraulic
Lifting capacity:	Max. 200kg
Mains voltage:	230 VAC, 50/60 Hz 120 VAC, 60 Hz 100 VAC, 50/60 Hz
Power consumption:	Max. 2.0 A
Insulation:	Class II
Device Classification:	I
Applied parts:	Type BF
Humidity proof:	IP21

Motor force:	6 kN
Motor usage:	Max. 2 minutes on / 18 minutes off
Hydraulic pump:	10kN
Frame material:	Steel frames
Frame finish:	Powder coated white RAL 9016
Conformity:	IEC/EN 60601-1 IEC/EN 60601-1-2 EC directives 2017/745 (MDR) 2012/19/EC 2011/65/EU



Warranty:	Table 3 years, Classique upholstery, 1 year
-----------	--

## CUSHION UPHOLSTERY

Classique:	PVC, phthalate free, biocompatible, abrasion resistant and flame retardant.
------------	---

## CUSHIONS FINISH

Linear:	Single padded, 40 kg/m <sup>3</sup>
---------	-------------------------------------

## OPTIONS

Wheels	
Safety switch for foot- or hand control	

Ergomax face opening

## DIMENSIONS / FUNCTIONALITY

