

SPECIFICATIONS

Dimensions :	205 x 67 cm	Motor force:	6 kN
Weight :	±79 kg	Motor usage:	Max. 2 minutes on / 18 minutes off
Height adjustment :	50 – 100 cm electric	Hydraulic pump:	10kN
	52 – 99 cm hydraulic	Frame material:	Steel frames
Lifting capacity:	Max. 200kg	Frame finish:	Powder coated white RAL 9016
Mains voltage:	230 VAC, 50/60 Hz	Conformity:	IEC/EN 60601-1
	120 VAC, 60 Hz		IEC/EN 60601-1-2
	100 VAC, 50/60 Hz		EC directives 2017/745 (MDR)
Power consumption:	Max. 2.0 A		2012/19/EC
Insulation:	Class II		2011/65/EU
Device Classification:	I	Warranty:	Table 3 years,
Applied parts:	Type BF		Classique 2025 upholstery,
Humidity proof:	IP21		2 years

CUSHION UPHOLSTERY

Classique:	PVC, phthalate free, biocompatible, abrasion resistant and flame retardant.
------------	---

CUSHIONS FINISH

Linear:	Single padded, 40 kg/m³
---------	-------------------------

OPTIONS

Wheels	Face opening
Safety switch for foot- or hand control	Ergomax face opening

DIMENSIONS / FUNCTIONALITY

