

## SPECIFICATIONS

Dimensions:	205 x 120 cm	Motor force:	8 kN
Weight:	±114 kg	Motor usage:	Max. 2 minutes on/ 18 minutes off
Height adjustment:	± 49 - 99 cm electric	Frame material:	Steel frames
Lifting capacity:	Max. 200 kg	Frame finish:	Powder coated white RAL 9016
Mains voltage:	100 - 240 VAC, 50/60 Hz	Conformity:	IEC/EN 60601-1 IEC/EN 60601-1-2 EC directives 2017/745 (MDR) 2012/19/EC 2011/65/EU
Power consumption:	Max. 3.0 A	<b>CE</b>	Warranty:
Insulation:	Class II		Table 3 years, Classique upholstery, 1 year
Device Classification:	I		
Applied parts:	Type BF		
Humidity proof:	IP21		

## CUSHION UPHOLSTERY

Classique: PVC, phthalate free, biocompatible, abrasion resistant and flame retardant.

## CUSHION FINISH

Linear: Single padded, 40 kg/m<sup>3</sup>

## OPTIONS

Wheels  
Cushion width 100cm

Face opening  
Ergomax face opening

## DIMENSIONS / FUNCTIONALITY

