

SPECIFICATIONS

Dimensions:	205 x 120 cm	Motor force:	8 kN
Weight:	±114 kg	Motor usage:	Max. 2 minutes on/ 18 minutes off
Height adjustment:	± 49 - 99 cm electric	Frame material:	Steel frames
Lifting capacity:	Max. 200 kg	Frame finish:	Powder coated white RAL 9016
Mains voltage:	100 - 240 VAC, 50/60 Hz	Conformity:	IEC/EN 60601-1 IEC/EN 60601-1-2 EC directives 2017/745 (MDR) 2012/19/EC 2011/65/EU
Power consumption:	Max. 3.0 A	CE	Warranty: Table 3 years, Classique 2025 upholstery, 2 years
Insulation:	Class II		
Device Classification:	I		
Applied parts:	Type BF		
Humidity proof:	IP21		

CUSHION UPHOLSTERY

Classique:	PVC, phthalate free, biocompatible, abrasion resistant and flame retardant.
------------	---

CUSHION FINISH

Linear:	Single padded, 40 kg/m³
---------	-------------------------

OPTIONS

Wheels	Face opening
Cushion width 100cm	Ergomax face opening

DIMENSIONS / FUNCTIONALITY

