


## SPECIFICATIONS

|                            |                           |   |  |
|----------------------------|---------------------------|---|--|
| Dimensions:                | 206 x 67 cm               | Motor force:  | 8 kN   |
| Weight:                    | 99 kg                     | Motor usage:  | Max. 2 minutes/ 18 minutes off   |
| Height adjustment:         | ± 47 – 97 cm              | Hydraulic pump:   | 10 kN  |
| Lifting capacity:          | Max. 200 kg               | Frame material:   | Steel frames   |
| Mains voltage:             | 100 - 240 VAC, 50 - 60 Hz | Frame finish:   | Powder coated white RAL 9016   |
| Current:                   | Max. 3.0 A                | Conformity:   | IEC/EN 60601-1   |
| Standby power consumption: | Max. 3 W                  |  | IEC/EN 60601-1-2   |
| Insulation:                | Class II                  |   | EU Regulation 2017/745   |
| Device Classification      | I                         |   | 2012/19/EU WEEE  |
| Applied parts              | Type BF                   |   | 2011/65/EU RoHS  |
| Humidity proof:            | IP21                      | Warranty:   | Table 3 years,<br>Premium upholstery 3 years,<br>Classique upholstery 1 year |

Note: motor adjustment is slower when running on battery.

## CUSHIONS UPHOLSTERY

|            |   |
|------------|---|
| Classique: | PVC, phthalate free, biocompatible, abrasion resistant and flame retardant. |
| Premium:   | PU, phthalate free, biocompatible, abrasion resistant and flame retardant.  |




## CUSHIONS FINISH

|              |   |
|--------------|---|
| Linear:      | Single padded, 40 kg/m <sup>3</sup>                         |
| Progressive: | Dual padded, 100 kg/m <sup>3</sup> and 40 kg/m <sup>3</sup> |

## OPTIONS

|                   |   |                        |   |
|-------------------|---|------------------------|---|
| Ergomax:          | PU foam, closed skin, low density<br>5° Shore   | Memory function:       | Three programmable settings for height adjustment<br>Bluetooth wireless operation         |
| Cushion heating:  | 34 °C, 36 °C and 38 °C<br>Bluetooth wireless operation                                | LED Pack:              | Chassis illumination 14 colours, including brightness control<br>Power consumption = 10 W |
| Battery function: | Automatic timer 30 minutes<br>Maximum consumption = 180 W<br>Li-ion 25,9 VDC, 2,25 Ah | Safe Motion Detection: | Power consumption = 1 W   |

FUNCTIONALITY/ DIMENSIONS



Model also available without flexion T5X

