

# Product datasheet

## G2 Trio

*The all-rounder of the Gymna G2 series, the G2 Trio treatment and examination couch has an adjustable head- and back section, making it the perfect solution for treatments in sitting and lying position. It is truly versatile due to its simplicity. The G2 Trio can be equipped with armrests for optimal relaxation in the prone position.*

### Standard features

- Three-section couch with adjustable head-, and back section
- All sections supported by gas springs
- Electric height adjustment with foot switch
- Cushions with 40 mm padding
- 15 cm clearance underneath frame for use with patient hoist (not compatible with 150 007)
- Non slip feet
- Cushions in a choice of 11 colours

### Ordering data

160 040 G2 Trio

### Optional


- 150 029 G2 Hand control
- 150 007 G2 Peripheral switch
- 150 001 G2 Hydraulic
- 150 062 G2 Wheel mechanism
- 150 002 G2 Arm supports\*
- 150 003 G2 Cushion with face opening
- 150 004 G2 Cushions in non-standard colour

### Accessories

- 150 106 G2 Paper roll holder 50 cm
- 150 117 G2 Paper roll holder 60 cm
- 337 018 Paper rolls 50 cm (per 12 pieces)
- 332 915 Foot end protector for cushion
- xxx xxx Terry-cloth cover, stretch material (available in 9 colours 'Style')

\* Automatically includes 150 003

### Specifications

Length x width:	± 205 x 67cm
Height adjustment:	± 47 – 100 cm
Weight:	± 80 kg
Max. lifting capacity:	200 kg
Electric motor:	6.000 N
Hydraulic pump:	10.000 N
Power supply:	230 VAC, 50-60 Hz, 1,0 A 120 VAC, 60 Hz, 1,7 A
Insulation class:	II
Moister resistance:	IP 21
Classification:	I
Conformity:	IEC/EN 60601-1 IEC/EN 60601-2 MDD 93/42/EEC
	
Guarantee:	3 years (upholstery 1 year)



160 040



**G2 Trio**



G2 Trio shown with options:

- Peripheral switch
- Wheel mechanism
- Arm supports

**Functionality/dimensions**



**Cushion colours**

**Style**



**Fashion**

