

gymna[®]

SPECIAL EDITION

gymna.GO Carbon Black



Basically Perfect.
Now Boldly Black.



Discover the gymna.GO

Carbon Black Edition

Gymna's most popular treatment table is now available in an exclusive Carbon Black Edition.

The robust black frame and matching upholstery give your practice a powerful, modern look without compromising on comfort and ease of use.

Select your version:



gymna.GO D1



gymna.GO T3x



gymna.GO T5x

Design that stands out.
Performance that convinces.



Carbon Black



Why choose the **Carbon Black Edition?**

- **Design with character:** the dark finish exudes professionalism and style.
- **Easy to move:** four swivel castors with locks make moving your treatment table smooth and safe.
- **Smart height adjustment:** With the innovative i-Control system, you can adjust the height around the table without interrupting your work.
- **Hygienic comfort:** The Ergomax face cushion is waterproof, washable and ergonomically shaped for maximum comfort.

Comfort you feel. Quality you notice.

Made for intensive use

The gymna.GO Carbon Black Edition is manufactured in the award-winning “factory of the future”. With robot-welded frames, durable powder coating and maintenance-free bearings, you get a table that will continue to operate smoothly and quietly for years to come, even in busy practices.

Materials that make a difference

The cushions are covered with premium faux leather in a deep black colour scheme with a leather look and feel. The durable PE foam filling provides optimal support and comfort, even during prolonged treatments.





Made to
perform every day



Standard included **Ergomax face cushion**

Maximum comfort, optimal hygiene

The Ergomax face cushion offers maximum hygiene and comfort:

- Heart-shaped opening for free breathing
- Soft, seamless finish
- Double-sided use, fits every body type
- Water-resistant, removable and easy to clean






Ready to GO?
Now boldly black.


Specifications

Feature	gymna.GO Carbon Black
Dimensions	205 x 67 cm (l x w)
Weight	circa 78 kg (D1) / 79 kg (T3x) / 81 kg (T5x)
Height adjustment	50 – 100 cm
Lifting capacity	Max. 200 kg
Motor power	6 kN
Motor use	2 min on / 18 min off
Mains voltage	230 / 120 / 100 VAC
Power consumption	Max. 2,0 A
Frame	Steel, powder coating RAL 9004 (black)
Cushion finish	Linear, single padded 40 kg/m ³
Upholstery	Premium PU Ink black
Device classification	Class I
Warranty	3 years on the table, 3 years on the upholstery


gymna.GO D1



gymna.GO T3x



gymna.GO T5x



gymna.GO

Carbon Black Edition



Need advice? Contact a Gymna dealer near you or visit gymna.com.

GymnaUniphy NV Pasweg 6A | B-3740 Bilzen, Belgium

Tel: +32 (0)89 510 532



Want to know more or
order directly?
Scan the QR code

www.gymna.com

gymna 12°