

Expand your imagination

ONE



gymna[®]

Expand your imagination

Our imagination shapes everything we do.

It helps us to express ourselves, unveiling who we really are. Both your practice as your treatment table should reflect the real you. Venture into the imaginary; expand your imagination even further. Immerse yourself into the gymna.ONE universe; dare to dream big and awake that imagination in your practice.





Refined design

gymna.ONE exudes design, the colours, materials and textures resonate nature. The use of durable and natural materials, combined with proven Gymna quality and reliability make it unparalleled. gymna.ONE emanates the harmonious atmosphere you imagined, it truly is the treatment table you aspire to have in your practice.



Your companion in daily practice

gymna.ONE is your companion in daily practice, its refined design combined with first-rate features optimize conditions for your treatments. Features such as i-Control, the intelligent foot-operated height adjustment system, the innovative Ergomax face support, positionable in two directions, and the lateral supports, that empower you to work closer to the patient. Choosing gymna.ONE as your trusted work companion, will benefit both you and your patients.



Ergomax face support

Made out of soft, seamless and comfortable foam, it is ergonomically shaped. The Ergomax face support is removable, waterproof and washable, allowing each patient to lie on a table that feels fresh and hygienic.



i-Control

The intelligent foot-operated height adjustment system for easy, accurate and safe height adjustment. Offering two operation modes, for continuous and pulsed adjustments.



Lateral supports

Allowing you even more freedom of movement while maintaining the ideal, ergonomically correct working position. They facilitate working close to the patient, without causing stress on shoulders or back.

Trust and confidence for your patient

Your patients are your prime focus. You want to offer an incomparable and foremost soothing and pleasant experience. Your goal is having your patients leave your practice feeling content and positive. The gymna.ONE's comfort and style will assist you in successfully building trust and erasing any of your patients' anxieties. Patients will undergo their next treatment session calmly and fully relaxed. As a result, you will have access to optimal conditions for a successful treatment.

Every gymna.ONE design detail has been developed with one goal in mind: the patient's comfort.



Remo

gymna.C
cushion
the cushi
level suite



Remote cushion heating

gymna.ONE is equipped with our remote cushion heating. The heating pads within the cushions generate a comfortable warmth level suited to your patient's individual needs.



Arm supports

Long arm supports provide ample space to let the arms rest completely. Making your patients feel secure and totally relaxed.



Body-shaped cushions

The body-shaped cushions add extra stability for your patients. Offering them a pleasant and comfortable treatment experience.



Ink Black



Ash Black



Midnight Blue



Warm taupe



Ruby Red



Caramel Latte



Ivory White



Dark Chocolate

Colours

Personalise the design of your gymna.ONE. Choose one of the 8 colours that reflect your professionalism and create the right atmosphere for your practice.





gymna.ONE *Carbon Black Edition*

Luxuriate even more in the beautiful design of the gymna.ONE with the special Carbon Black edition. In this timeless and classic finish, the gymna.ONE exudes style and refinement, the colours, materials and textures resonate nature. The use of durable and natural materials, combined with proven Gymna quality and reliability make it unparalleled. gymna.ONE emanates the stylish atmosphere you imagined, this is the treatment table that truly stands out in your practice.

Pure design

The classic carbon black colour of the gymna.ONE Carbon Black edition's cushions, alongside the beautiful shape of the aluminium lifting arms and the black coloured closed base make this gymna.ONE version exceptionally unique and memorable. The closed base in high quality beech plywood is finished with satin black laminate. The curved and flat surfaces of the base reflect and coruscate the surrounding light, creating a dynamic visual effect while highlighting the different shades of black.



Specifications

Dimensions	± 205 x 67 cm
Height adjustment	± 50 – 100 cm
Weight	± 127 kg
Max. lifting capacity	Max. 200 kg
Electric motor	8,000 and 6,000 N
Cushion temperature	34 – 36 – 38°C
Power supply	100 – 240 VAC, 50/60 Hz
Current consumption	2.3 A
Usage	Max. 2 minutes on / 18 minutes off
Insulation class	II
Moisture resistance	IP21
Classification	I
Guarantee	3 years (upholstery 3 years)
Conformity	IEC/EN 60601-1 IEC/EN 60601-1-2 EC directives 2017/745 (MDR) 2012/19/EC 2011/65/EU



Configuration




gymna.ONE configuration	T7	Q8
Features		
Wheels	●	●
Safety switch	●	●
i-Control	●	●
Gymna Stability Profile	●	●
Choice between black or white	●	●
Top Sections		
Ergomax face opening	●	●
Arm supports	●	●
Lateral supports	●	●
Manual lumbar flexion	●	-
Motorized lumbar flexion	○	-
Motorized back section	-	●
Motorized sitting section	-	●
Cushions		
Progressive cushions	●	●
Remote cushion heating	●	●
Premium upholstery	●	●
Width 67 cm	●	●

● Standard ○ Option - Not available



Your dealer



gymna.com

gymna[®]

GymnaUniphy NV Pasweg 6A | B-3740 Bilzen, Belgium

Tel: +32 (0)89 510 532